CO₂ Chamber for mice and rats

Simple, Quick & Easy to Use!

Door

- Place animals in chamber without losing CO2.
- Provides for ultra quick animal

The Cleanest CO2 Induction Chamber on the Market!

Delivery Tube Within Chamber

Fills chamber from bottom, decrease induction time, "more humane".



Quick Connect w/valve

Closes when disconnected, keeps CO₂ in tube and chamber.

7' Tube

Easily connects to CO₂ source.

Clear Acrylic Top

Easy to see animals in chamber.



Solid Sides

Prevents viewing of animals during the process.

Small Foot Print

9" diameter, easy to keep on bench.





Disposable Bag Liner

- Animal feces & urine remain in bag with animal.
- Chamber is always clean.

CO₂ Chamber Cat #: COIC

Includes: See thru acrylic lid with door & quick connect w/ valve, 7' tubing, 50 disposable insert bags.



Cat #: COBL Disposable Bag Liner (50 bags/pk)

Use w/ CO2 Chamber.



Animal Identification and Marking Systems, Inc.

Voice 607.324.6752 • www.animalid.com

All text and graphic are © 2012 by Animal Identification and Marking Systems, Inc., All rights reserved. No text or graphics may be copied or used without the written permission from Animal Identification and Marking Systems, Inc.