

"Single Temp Warming/Induction Chamber"




Post-Surgery Warming Chamber * Pre-Surgery Warming Induction Chamber

Warming/Induction Chamber provides continuous warmth to animals being anaesthetized before surgery or recovering from surgery.




Filter Top: Allows animal to fully recover in chamber while unattended. Includes membrane for sealing filter top and outlet port to exhaust gas when used as an induction chamber.




Location of Warming Elements Under Stainless Steel Plate

Small Chamber Full-Floor Warming




Cat #
WCS0 / WCS4
WCS8 / WCS9

Large Chamber, Full-Floor Warming




Cat #
WCL0 / WCL4
WCL8 / WCL9

Large Chamber, Half-Floor Warming



Cat #
WCH0 / WCH4
WCH8 / WCH9

Cool side of cage.



Cat #
WCS5

Small Chamber Half-Floor Warming

Important Note: Cage floor temperatures were measured with an ambient room temperature of 70°F (21°C). Room conditions and table temperature will slightly increase/decrease cage floor temperature.

Warming Chamber

Std. Small Shoe Box Cage, 7" x 11" x 5", Std. Large Shoe Box Cage, 10" x 19" x 8".
Inlet/Outlet Port for Gas Anesthetic When Used As An Induction Chamber.
Heavy Duty Power Packs Provide Preset Temp. Universal Input (110-220 V) for use in any country.
Warming Element Sealed Under Stainless Steel Floor. Bedding Material Can Be Used In Cage.

Simple and Easy to Use Plug-in Power Supply



Cat#

Description

Small Chamber 7" x 11" x 5", Full-Floor Warming

WCS0 ~30°C Single Temp Small Warming/Induction Chamber
WCS4 ~34°C Single Temp Small Warming/Induction Chamber
WCS8 ~38°C Single Temp Small Warming/Induction Chamber
WCS9 Variable Temp Small Warming/Induction Chamber*

Small Chamber 7" x 11" x 5", Half-Floor Warming

WCS5 Variable Temp Small Warming/Induction Chamber*

Large Chamber 10" x 19" x 8", Full-Floor Warming

WCL0 ~30°C Single Temp Large Warming/Induction Chamber
WCL4 ~34°C Single Temp Large Warming/Induction Chamber
WCL8 ~38°C Single Temp Large Warming/Induction Chamber

Large Chamber 10" x 19" x 8", Half-Floor Warming

WCH0 ~30°C Single Temp Large Warming/Induction Chamber (H)
WCH4 ~34°C Single Temp Large Warming/Induction Chamber (H)
WCH8 ~38°C Single Temp Large Warming/Induction Chamber (H)

Chamber Include: UIHD-Power Pack, Induction Conversion Kit & Operating Instructions
*~30°C, ~36°C, ~38°C, ~42°C, ~49°C
(actual temperatures will be dependent on your room conditions)

Extra Power Packs

Universal Input 110-220V, can be used with any warming chamber

PP30 ~30°C Single Temp, UIHD-Power Pack
PP34 ~34°C Single Temp, UIHD-Power Pack
PP38 ~38°C Single Temp, UIHD-Power Pack
PPVT Variable Temp, UIHD-Power Pack (30°C, 36°C, 38°C, 42°C, 49°C)