

**New**

# CO<sub>2</sub> Chamber

for mice and rats

## Door

- Place animals in chamber without losing CO<sub>2</sub>.
- Provides for ultra quick animal induction.



***Simple, Quick & Easy to Use!***

*The Cleanest CO<sub>2</sub> Induction Chamber on the Market!*

## Delivery Tube Within Chamber

- Fills chamber from bottom, decrease induction time, "more humane".

## Quick Connect w/valve

- Closes when disconnected, keeps CO<sub>2</sub> in tube and chamber.

## 7' Tube

- Easily connects to CO<sub>2</sub> source.

## Clear Acrylic Top

- Easy to see animals in chamber.



## Solid Sides

- Prevents viewing of animals during the process.

## Small Foot Print

- 9" diameter, easy to keep on bench.



## Disposable Bag Liner

- Animal feces & urine remain in bag with animal.
- Chamber is always clean.

## Cat #: COIC CO<sub>2</sub> Chamber

Includes: See thru acrylic lid with door & quick connect w/ valve, 7' tubing, 50 disposable insert bags.



## Cat #: COBL Disposable Bag Liner (50 bags/pk)

Use w/ CO<sub>2</sub> Chamber.



For customer service, call 1-800-766-7000.  
To fax an order, use 1-800-926-1166.  
To order online: [www.fishersci.com](http://www.fishersci.com)