



# Dissection Platform

**TWO  
Ultra-Bright  
LED Lights**



**Magnifying  
Glass**  
**main lens: 2.3X**  
**inset lens: 4.0X**




Cat#: DTW1

AIMS Blue Wax Dissection Tray  
(9" x 6.5" x 1")

**USE:** Can be used for dissections, surgery, micro manipulations, general micro-utility platform (i.e. organ inspections, dissections, etc.).

- Two Flexible LED lights - provides direct illumination with no shadows and no heat.
- Flex arm, non-obstructive rimless magnifying glass (main lens: 2.3X, inset lens: 4.0X). (Allows continuous eye contact with object when moving into magnification area).
- Limb restraints.
- Non-skid feet for added stability.
- Universal input power supply (100 – 240 VAC, 50-60 Hz).
- Size: 9" X 13.5"

**ORDER: ASS7**



## Tissue Processing Disc

Process tissue using a clean blue wax surface.

« layout, pin, & trim tissue »

Use for processing, inspecting, & documenting tissue pathology.

Cat#: LTPD  
Cat#: MTPD  
Cat#: STPD  
Cat#: TPIN

**Large Tissue Processing Disc** (4" dia., ¼" thick, 10/pk)  
**Medium Tissue Processing Disc** (2 ¾" dia., ¼" thick, 10/pk)  
**Small Tissue Processing Disc** (2" dia., ¼" thick, 10/pk)  
**Ball Head Tissue Pins, ss** (1½" long, 50/pk)



Animal Identification and Marking Systems, Inc., V: 607.324.6752 • www.animalid.com • aims@animalid.com • F: 607.324.6753

All text and graphic are © 2012 by Animal Identification and Marking Systems, Inc., All rights reserved. No text or graphics may be copied or used without the written permission from Animal Identification and Marking Systems, Inc.