
Surgery Platforms

Surgery Platform



USE: Can be used for surgery, micro manipulations, general micro-utility platform (i.e. organ inspections, dissections, etc.).


- <u>Two</u> Flexible LED lights provides direct illumination with no shadows and no heat.
- Flex arm, non-obstructive rimless magnifying glass (main lens: 2.3X, inset lens: 4.0X). (Allows continuous eye contact with object when moving into magnification area).
- Limb restraints.
- Non-skid feet for added stability.
- Universal input power supply (100 240 VAC, 50-60 Hz).
- Size: 9" X 13.5"

ORDER: <u>ASS7</u>

Thermo Recovery Pad

Animal Recovery Made Simple and Easy !





Ultra-Safe Thermo Regulation:

- will not over heat ani-
- will not burn animal
- no electricity used

Place animal on pad for recovery. Recovery pad will provide over 2 hrs. of Thermo-Regulation.

Thermo Surgery Platform



- <u>Two</u> Flexible LED lights provides direct illumination with no shadows and no heat. (now up to 50% brighter)
- Flex arm, non-obstructive rimless magnifying glass (main lens: 2.3 X, inset lens: 4.0 X). (Allows continuous eye contact with object when moving into magnification area)
- Limb restraints.
- 4 Space Gels for continuous heat.
- 2 Space Blankets to retain heat.
- Thermo Reflector and Turbo Heat Sink to maximize heat transfer.
- Non-skid feet for added stability.
- Universal input power supply.
- Size: 13.5" (L) x 9" (W)

ORDER: ASS7T

ORDER: TCRP

- Convenient size 13.5" (L) x 7.5" (w)
- 3 Space Gels for continuous heat.
- 4 Space Blankets to retain heat.
- Thermo Reflector and Turbo Heat Sink to maximize heat transfer.
- Non-skid feet for added stability.
- Fits into a standard 18 (L) x 9 (W) shoe box cage.

